



Solid Rock

José R. Valverde
Scientific Computing
CNB/CSIC

jrvalverde@cnb.csic.es
CC-BY-NC-SA

Making Movies: Chimera

- Basic commands
- Spin/rock
- Center
- Increase quality
- Creating animations
- Volume rendering
- Morphing

Making it simple

- Open chimera
 - Load molecule(s)
 - Tools → Utilities → Movie recorder
 - Record
 - Move molecule
 - Stop
 - Make movie



?

Automatic movement

- Favorites → Command Line
 - **turn y 3 120**
- Tools → Utilities → Movie recorder
 - Record
 - Press enter on command line
 - When done press Stop and Make Movie



?

Basic commands

- You can use several commands separated by a ';'
 - movie record
 - movie stop
 - movie encode
- General form: **movie** *action options*
 - **movie record directory** *dir1/dir2/*
supersample 3 size w,h raytrace true
 - **movie encode mformat** *mp4 output mymv*

Automatic recording

- The problem lies in hand switching
- Let us use the command line:
 - movie record ; turn y 3 120 ; wait; movie stop
 - movie encode output mymovie



?

Controlling movement


- You can use any visualization command between **movie** commands
- Remember to **wait** for a long command to finish before issuing the next
- Try rocking instead of spinning
 - movie record; rock y 3 68; wait; rock x 3 68; wait; movie stop; movie encode output mymovie.mov



?

Control animation

- You can use any visualization command
 - On the menu
 - On the command line
- Better if you define layout first
 - Change window size
 - Presets → Interactive 1
 - Select → residue
 - Actions → set pivot...



Control animation (2)

- `window size 640 480`
- `preset apply interactive 1`
- `cofr :GLN`
- `movie record ; rock y 3 68 ; wait ; rock x 3 68; wait ; movie stop ; movie encode output mymovie.mov`
- You can put everything in a script file **.cmd** and load it with **File → open**



?

Improving quality

- Record with supersampling (e.g. 3): draw at a bigger size and resize down
 - movie record supersample 3
- Record at desired size
 - movie record size 640,480
- Increase output bit rate for a given window size, 640x480: 200 (low), 1000 (medium), 6000 (high)
 - movie encode bitrate 10000
- Use appropriate encoding (e.g. mov, mp4)
- Add fill frames



Raytracing

- Very slow, may hang GUI for $> \sim 200$ frames
- Produces a POVray file for each frame and renders it using POVray

```
window size 640 480
```

```
~ribbon
```

```
show :GLN
```

```
repr sphere
```

```
color byelement
```

```
reset
```

```
set bg_color black
```

```
movie record supersample 3 raytrace true
```

```
wait 10 ;rock y 3 68; wait; rock x 3 68; wait; wait 10
```

```
movie stop
```

```
movie encode mformat mp4 bitrate 10000 output 08-raytrace.mp4
```




Animation

- Start by marking all needed positions:
 - savepos start
 - center :GLN; savepos centered
 - ~select 0; move x 30; move y 40; turn x 90; select 0; savepos unbound
- Chimera moves among positions with **reset**. If more than one frame is requested, linear interpolation will be used.

Prepare world

- Select initial position
- Select conformation/display
- Add labels:
 - 2dlabel create title “XX” xpos 0.3 ypos 0.9 color white
- Select background color
 - set bg_color gray
- etc...

You can do anything...

- While recording:
 - Add/remove surfaces
 - Fade in/out objects
 - Identify Hydrogen bonds
 - Color by any criteria
 - Change display modes (CPK, surface..)
 - Anything!

Getting it done

- Work on a script
 - Prepare positions
 - Set up initial scene
 - Comment out (#) movie record
 - Keep on adding commands until OK
 - Add movie stop; movie encode
 - Uncomment movie record
- Load script



Volume rendering

- You can also animate volumetric data
 - e.g. from EMSD:1321
 - open #0 emdbID:1321
 - volume #0 level 127
 - volume #0 color yellow
 - movie record
 - rock x 3 68; rock y 3 68
 - roll 1,1,0 3 68
 - movie stop; movie encode



Morphing

- You can morph among different structures too:
 - morph start #0.1
 - morph interpolate #0.2
 - ...
 - morph interpolate #0.10
 - morph movie
 - movie record
 - coordset #1 1, ; wait; movie stop; ...



Morphing example

- Fetch by ID: PDB:1az6
- Tools → Structure comparison → Morph conformations
- Add → Add all variants, then add again variant #1
- Action on create → Hide conformations
- Create
- Play







That's all, folks!

QUESTIONS?